# MICRO REPORT Open Access



# Axonal membrane stretch suppresses neuronal excitability by activating mechanosensitive K2P channels at the node of Ranvier

Hirosato Kanda<sup>1,2\*</sup>, Koichi Noguchi<sup>2</sup> and Yi Dai<sup>1,2,3</sup>

## **Abstract**

Saltatory conduction is the propagation of action potentials along myelinated nerves, which enables fast propagation through the node of Ranvier. Recently, we demonstrated that K2P channels, TWIK-related K+ channel-1 (TREK-1), and TWIK-related arachidonic acid-activated K+channel (TRAAK), are highly expressed in the mammalian node of Ranvier of sensory nerves and have an important role in action potential repolarization instead of voltage-gated K<sup>+</sup> channels. TREK-1/TRAAK channels are activated by membrane depolarization as well as various stimuli, such as temperature, pH, arachidonic acid, and mechanical membrane stretch. Although membrane mechanical stretch has been suggested to modulate action potential conduction, how membrane stretching modulates intrinsic electrophysiological properties at the node of Ranvier remains unclear. In the present study, we examined the effects of membrane stretch on neuronal membranes at the node of Ranvier in rat sciatic nerves. The single-channel conductance was approximately 90 pS at 80 mV. Membrane stretch increased the single-channel event numbers and open probability in a pressuredependent manner. Consistent with single-channel activity, intra-pipette positive pressure increased outward leak currents and decreased membrane excitability in a whole-cell configuration. Furthermore, blockage of TREK-1/TRAAK channels by Ba<sup>2+</sup> reversed the changes in the intrinsic electrophysiological properties induced by intra-pipette pressure. These results indicate that the activation of mechanosensitive TREK-1/TRAAK channels may suppress neuronal excitability following axonal stretch. Our findings suggest that TREK-1/TRAAK channels may play an important role in the prevention of ectopic action potential discharge at the axon by intense mechanical nerve stretch under physiological conditions.

**Keywords** Node of Ranvier, Two-pore-domain potassium channels, Membrane stretch, Action potential, Single-channel activity

Hirosato Kanda

kanda@hyo-med.ac.jp

### Introduction

The axon ensures stable action potential (AP) propagation from the axon initial segment to the axonal terminus [1]. While axons are generally considered reliable transmission cables, an AP can be initiated at the axon in an unusual condition in which intense pressure is directly applied to axons, or because of axonal injury, such as nerve compression syndrome and neuropathic pain [2, 3]. Recently, we and other have demonstrated



© The Author(s) 2023. **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit http://creativecommons.org/licenses/by/4.0/. The Creative Commons Public Domain Dedication waiver (http://creativecommons.org/publicdomain/zero/1.0/) applies to the data made available in this article, unless otherwise stated in a credit line to the data.

<sup>\*</sup>Correspondence:

School of Pharmacy, Hyogo Medical University, Kobe, Hyogo 650-8530,

 $<sup>^{2}</sup>$  Laboratory of Basic Pain Research, Hyogo Medical University, Kobe, Hyogo 650-8530, Japan

<sup>&</sup>lt;sup>3</sup> Department of Anatomy and Neuroscience, Hyogo Medical University, Nishinomiya, Hyogo 663-8501, Japan

Kanda et al. Molecular Brain (2023) 16:8 Page 2 of 4

that two tandem pore domain potassium (K+) channel (K2P channel) proteins, TREK-1 (TWIK-related K<sup>+</sup> channel), and TRAAK (TWIK-related arachidonic acidactivated K<sup>+</sup> channel) are highly expressed in the mammalian nodes of Ranvier (NRs) of the sensory nerve, and that AP repolarization at NRs relies on TREK-1/TRAAK channel activation [4, 5]. Interestingly, TREK-1/TRAAK are activated by membrane depolarization and by various stimuli, such as temperature, pH, arachidonic acid, and mechanical membrane stretch [5–8]. The suppression of TREK-1/TRAAK channel activities by cold temperature or pH delays AP repolarization and conduction velocity, and impairs the regeneration of high-frequency AP trains [4, 9]. Since TREK-1 and TRAAK are mechanosensitive K<sup>+</sup> channels, axonal membrane stretching is likely to modulate the intrinsic electrophysiological properties of myelinated nerves. Liu et al. recently developed an experimental technique to apply an axonal mechanical stretch to ex vivo sciatic nerve preparation. Using these methods, they demonstrated that mechanical nerve stretch instantaneously mediates conduction delay in A-fibers [10]. Although membrane mechanical stretch has been suggested to modulate AP conduction, it is still unknown how membrane stretching modulates the excitability of axonal membranes at the NRs. In this study, we performed patch-clamp recording on ex vivo sciatic nerve preparations and examined the effects of membrane stretch on single-channel activity and intrinsic electrophysiological properties of the NRs.

# **Results and discussion**

To determine the effect of membrane stretch at NRs of myelinated sensory nerves, we performed a pressureclamped patch-clamp on intact NRs of rat sciatic nerve [11] (Fig. 1A). We first analyzed single-channel activities at the NRs by applying various negative pressures on the axonal membrane with a high-speed pressure-clamp device under the cell-attached configuration (Fig. 1B). The single-channel conductance was approximately 90 pS at 80 mV, which was consistent with the channel conductance of the heteromeric TREK-1/TRAAK channels [4] (Fig. 1C). The single-channel conductance was not affected by membrane stretch (Fig. 1C). The single-channel event numbers and open probability increased in a pressure-dependent manner and showed a significant difference from the negative pressure of 60 mmHg or more (Fig. 1D, E). These results show that TREK-1/TRAAK channels are highly sensitive to mechanical membrane stretch, and suggest that TREK-1/TRAAK activation by membrane stretch may modulate neuronal excitability at


At NRs, strong leak potassium outward (lK<sub>leak</sub>) currents following depolarizing voltage steps mainly mediated by

activation of TREK-1/TRAAK channels [4]. To determine whether membrane stretch enhances the IK<sub>leak</sub> currents by this activation at the NRs, we applied positive pressure into cells through a recording pipette under the whole-cell configuration (Fig. 1F). By applying positive pressure in this configuration, we applied membrane tension in the same direction from inside to outside of the cell as cell-attached configuration. Because the seal between the cell membrane and patch pipette is easily broken by internal positive pressure, we applied a maximum pressure of 10 mmHg for the experiment. This positive pressure enhanced the outward currents in response to the voltage steps (Fig. 1F). Consistent with singlechannel activities, intra-pipette positive pressure significantly increased IK<sub>leak</sub> conductance compared to that at 0 mmHg (Fig. 1G). To confirm whether the enhanced outward currents were mediated by TREK-1/TRAAK channels, we pharmacologically blocked these channels by applying Ba<sup>2+</sup>. The bath application of Ba<sup>2+</sup> significantly suppressed pressure-enhanced outward currents and IK<sub>leak</sub> conductance (Fig. 1G, H). These results clearly show that intra-pipette pressure enhances K+conductance by activating TREK-1/TRAAK channels at the NRs.

We further confirmed whether membrane stretch modulates nodal membrane excitability in NRs in the current-clamp mode. The input resistance was significantly decreased by 10 mmHg intra-pipette positive pressure (Fig. 1I, J). The resting membrane potential was significantly decreased (Fig. 1K) and the rheobase was significantly increased by an intra-pipette pressure of 10 mmHg (Fig. 1L). The blockage of TREK-1/TRAAK channels by Ba<sup>2+</sup> reversed the changes in the intrinsic electrophysiological properties induced by intra-pipette pressure (Fig. 1J–L).

In conclusion, our study demonstrates that mechanosensitive TREK-1/TRAAK channels may be activated by axonal stretch. The activation could suppress neuronal excitability. We previously found that TREK-1 and TRAAK are highly expressed in the NRs of mammalian sensory nerves, which leads to AP repolarization instead of voltage-gated K<sup>+</sup> channels [4]. In the present study, single-channel activities with channel conductance of 90 pS (TREK-1/TRAAK-like) were increased by membrane stretch in a pressure-dependent manner. Consistent with single-channel activity, intra-pipette positive pressure increased IK<sub>leak</sub> currents and suppressed membrane excitability. Furthermore, pharmacological blockage of TREK-1/TRAAK channels reversed the changes in intrinsic electrophysiological properties induced by intra-pipette pressure. These results indicate the importance of TREK-1/TRAAK channels in the prevention of ectopic AP discharge at the axon by intense mechanical nerve stretch under physiological conditions.

Kanda et al. Molecular Brain (2023) 16:8 Page 3 of 4



**Fig. 1** Effects of membrane stretch of neuronal membrane at the node of Ranvier (NR) of rat sensory nerve fiber. **A** Bright-field image showing patch-clamp recordings at an NR of Aβ-afferent nerve in ex vivo sciatic nerve preparation. **B** Left: recording setting of single-channel activity. Right: traces illustrating  $|K_{leak}|$  single-channel currents at different intensities of membrane stretch on the nodal membrane. Summary of  $|K_{leak}|$  single-channel conductance (**C**), single-channel events (**D**), and open probability (**E**) at different intensities of membrane stretch (n = 7). **F–H** Left: recording setting under the whole-cell patch-clamp. Right: Traces illustrating currents recorded at NRs by voltage steps (**F**), |-V| curves of inactivating currents (**G**), and whole-cell conductance of  $|K|_{leak}|$  (**H**) following intra-pipette application of 0 mmHg (n = 9) and 10 mmHg pressure (n = 5), 0 mmHg pressure with  $|B|^2$  (n = 5), and |A| 10 mmHg pressure with  $|B|^2$  (n = 5). Traces illustrate typical action potentials recorded at NRs (**I**), input resistance (**J**), resting membrane potential (**K**), rheobase (**L**) following intra-pipette application of 0 mmHg (n = 7) and |A| 10 mmHg pressure (n = 5), 0 mmHg pressure with  $|B|^2$  (n = 5), and |A| 10 mmHg pressure with  $|B|^2$  (n = 5). Arrow indicates NR. Data represent mean |A| 5EM; ns, not significantly different; |A| 20.05, |A| 20.01, |A| 3 mmHg pressure with Tukey's multiple comparisons test

Kanda et al. Molecular Brain (2023) 16:8 Page 4 of 4

As a limitation of this study, we used intra-pipette pressure to study ion channel function at the NRs instead of axonal stretch, because the membrane seal for patch-clamp recording is easily disrupted by membrane displacement. Thus, the application of mechanical pressure with a high-speed pressure-clamp device may not fully represent axonal stretch.

In conclusion, our study demonstrates the effect of membrane stretching on the intrinsic electrophysiological properties of NR. The findings provide important insights into the pathology of diseases, such as demyelination.

### Abbreviations

AP Action potential

K2P channel Two tandem pore domain potassium channel

 $\mathsf{IK}_{\mathsf{leak}} \qquad \qquad \mathsf{Leak} \ \mathsf{potassium} \ \mathsf{outward}$ 

NRs Nodes of Ranvier

TRAAK TWIK-related arachidonicacid-activated K<sup>+</sup> channel

TREK-1 TWIK-related K<sup>+</sup> channel

# **Supplementary Information**

The online version contains supplementary material available at https://doi.org/10.1186/s13041-023-01000-6.

Additional file 1. Materials and methods.

# Acknowledgements

We would like to thank Seiko Oki for technical support and Editage (www. editage.com) for the English language editing.

# **Author contributions**

Conceptualization: HK; performance of the experiments: HK; formal analysis: HK; Writing-original draft: HK; supervision: KN and YD. All authors read and approved the final manuscript.

### **Funding**

Not applicable.

# Availability of data and materials

All data generated during this study are available from corresponding author on reasonable request (Additional file 1).

# **Declarations**

### Ethics approval and consent to participate

All procedures were approved by Hyogo University of Health Sciences Committee on Animal Research. Consent to participate is not applicable for this study.

# Consent for publication

Not applicable.

### Competing interests

The authors declare no conflict of interest.

Received: 15 November 2022 Accepted: 10 January 2023 Published online: 17 January 2023

### References

- Debanne D, Campanac E, Bialowas A, Carlier E, Debanne D, Campanac E, et al. Axon. Physiol Physiol Rev. 2011;91:555–602.
- Chul Han H, Hyun Lee D, Mo Chung J. Characteristics of ectopic discharges in a rat neuropathic pain model. Pain. 2000;84:253–61.
- 3. Finnerup NB, Kuner R, Jensen TS. Neuropathic pain: Frommechanisms to treatment. Physiol Rev. 2021;101:259–301.
- Kanda H, Ling J, Tonomura S, Noguchi K, Matalon S, Gu JG. TREK-1 and TRAAK are principal K + channels at the nodes of Ranvier for Rapid Action potential conduction on mammalian myelinated afferent nerves. Neuron. 2019;45:67.
- Brohawn SG, Su Z, MacKinnon R. Mechanosensitivity is mediated directly by the lipid membrane in TRAAK and TREK1 K + channels. Proc Natl Acad Sci U S A. 2014;111:3614–9.
- Kang D, Choe C, Kim D. Thermosensitivity of the two-pore domain K + channels TREK-2 and TRAAK. J Physiol. 2005;564:103–16.
- Schewe M, Nematian-Ardestani É, Sun H, Musinszki M, Cordeiro S, Bucci G, et al. A non-canonical voltage-sensing mechanism controls gating in K2P K + channels. Cell. 2016;164:937–49.
- Levitz J, Royal P, Comoglio Y, Wdziekonski B, Schaub S, Clemens DM, et al. Heterodimerization within the TREK channel subfamily produces a diverse family of highly regulated potassium channels. Proc Natl Acad Sci U S A. 2016;113:4194–9.
- Kanda H, Tonomura S, Gu JG. Effects of cooling temperatures via thermal k2p channels on regeneration of high-frequency action potentials at nodes of ranvier of rat ab-afferent nerves. eNeuro. 2021;8:1–16.
- Liu J, Ganeshbabu N, Shalaby N, Chen L, Guo T, Feng B. Targeting twopore-domain potassium channels by mechanical Stretch instantaneously modulates action potential transmission in mouse sciatic nerves. ACS Chem Neurosci. 2021;12:3558–66.
- Kanda H, Tonomura S, Dai Y, Gu JG. Protocol for pressure-clamped patchclamp recording at the node of Ranvier of rat myelinated nerves. STAR Protoc. 2021;2:100266.

## **Publisher's Note**

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

# Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

### At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

